

45155-002



45155-202



45155-202
45155-002



45155-602



45155-702



AUF EINEN SICHEREN STAND DES GERÄTES ACHTEN.

*Ensure that the unit is in a secure position.
Veuillez à ce que l'appareil soit en position sûre.*

ÜBERPRÜFEN SIE REGELMÄßIG DEN FESTEN SITZ DER UNIVERSALKLEMME, GEGEBENENFALLS NACHZIEHEN.

*Regularly check that the multi clamp is firmly seated and secured, retighten if necessary.
Vérifiez régulièrement que la pince multi est bien en place et fixée, resserrez si nécessaire.*